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The interface includes a top menu bar with options like Camera, JoyStick, Sample Stage, Modes, Laser Power Adjust, Optical Simulator, JoyStick - Mask, Image Plane Mapping, Laser Spot Size Control, and I-V. Below this are several control panels: a 'Power Set to' panel with a digital display at 0 mW and a slider; a 'Connection Error - Pleas' panel; a 'Connection Error' panel with settings for Pol./Mech. Angle, Light Pol. Step Angle, Polarizer Steps, Analyzer Step Angle, Analyzer Steps, Wait, and Enable Acceleration; and an 'Objective' panel with a dropdown menu and Pmin/Pmax indicators.

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The MOKE & AHE interface features a 'Program Text' area on the left with a sequence of commands like 'Start', 'Set field to 10 mT', 'Set Hspeed to 400 mT/min', 'Stop measure', and 'End'. The main area contains several graphs: 'DC Graph', 'z-H Graph', 'B-H Graph', and 'R-H Graph'. A 'State Indicator' section shows 'Heater To Zero?' and 'Exciting Line' status. A 'Command Table' lists various actions like 'Set Temperature', 'Field', 'Current', 'Wait', 'Laser', 'Rotator #1', 'Rotator #2', 'Action', 'Hold', 'Clamp', 'Switch Heater', 'Measure', 'Variable Select', 'Stop Measure', 'Auto Functions', 'Focus', 'Calibrate B', and 'Calibrate c'. An 'Instruments Selector' panel on the right has buttons for Magnet, LIA DC, LIA f, LIA 2f, Read R, and Read I.

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The first graph shows the Curie temperature (Tc) vs. thickness for CrTe2, with regions for Paramagnetic (PM), Ferromagnetic (FM), and Ferromagnetic Insulator (FMI). The second graph shows Tc vs. thickness for Cr5Te8, with Paramagnetic and Ferromagnetic regions. The third graph shows the Kerr angle vs. magnetic field for FeTa3S6, highlighting 'Fresh charge' and '4 months aging' curves.

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The first graph shows Kerr rotation vs. temperature for a heterostructure, with curves for different magnetic fields (0 GPa, 10 GPa, 15 GPa). The second graph shows Kerr rotation vs. temperature for another heterostructure, with curves for different magnetic fields. The third graph shows Kerr rotation vs. thickness for Cr2Te3, with curves for different magnetic fields.

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